

biological · biophysical · engineering

moldaenke





Quick and reliable drinking water protection

the water watchdog

www.bbe-moldaenke.de



Automated online surveillance

The *Tox*Protect64 is an automated monitoring system for the protection of a drinking water supply against accidental or malicious contamination due to harmful substances. The characteristics of such threats are a relatively high concentration of dangerous substances that occurs suddenly. This challenging task requires a number of important criteria to be met. Such a system MUST be sensitive to a wide range of toxins. The detection of toxins using fish is a well-established method, with sensitivity data for nearly every combination of fish and toxin readily available.

The test organisms must react to substances harmful to humans. The use of fish gives the closest practical comparison available in the expected scenario. Falsealarms MUSTbe reduced to an absolute minimum in order to increase user-confidence and avoid unnecessary expense. The *Tox*Protect64 uses an integrated alarm verification system to fulfil this requirement. It is easyto operate at affordable costs since for a high level of security, it may be necessary to employ multiple monitoring locations.

The *Tox*Protect64 is a low-price, low-maintenance device which can be installed at various sites within a water company.



Operation

The *Tox*Protect64 user interface allows easycontrol of the instrument. All readings are shown as numbers or as charts. The parameters of the *Tox*Protect64 measuring processcan either be entered via the control panel or via the web interface.

Control panel

On the left, a parameter settings menu to select the relevant submenu, in the center a summary of parameter readings and on the right an LEDtest displaying the number of interruptions for each light barrier during the last 60s.



Detection Method

The ToxProtect64 monitors the swimming activity of up to 20 fish by measuring the frequency of interruptions of an array of light barriers. The result is given in interruptsper minute and fish. In the event of valuesfalling below a given threshold for a certain time, the alarm verification processis initiated.

Additionally, immobile fish at the bottom and in the upper region of the aquarium are registered.

The fish species employed is user-selectable, with recommendations given in the specifications. Best results are obtained using agile fish of 4-6 cm in length.

Alarm Verification

Due to natural random variations of fish behaviour, any alarm criterion may be reachedfrom time to time just by accident. Hence, to prevent false alarms, a verification system is required. This is achieved by increasing the illumination inside the aquarium during verification.

Normally, thisleadsto a dramatic increase in the activity of the fish. Under toxic conditions, this may not occur. The instrument thus accepts or rejects the alarm automatically.

Alarm Indication

F

Y

Ν

F

The alarm indicator (traffic light) shows the current operating mode:

Red:	alarm
Yellow:	alarm verification
Green:	normal operation
No light :	stand-by
Flashing green:	start-up
Flashing yellow:	failure



Instrument Malfunction Surveillance

The ToxProtect64 is equipped with internal sensorsto monitor and report on instrument malfunction including: Dinadequate sampleflow Ddrain blockage Dtemperature excursion Daccidental/unauthorised exposureof test chamberto ambient light Dhigh chlorine concentration Dloss of dechlorination reagent

Optional Dechlorination System

For chlorinated drinking water the ToxProtect64 can be equipped with a dechlorination unit with integrated level control, a peristaltic pump and a total chlorophyll electrode.



Verified Alarm of the ToxProtect64

Discarded Alarm of the ToxProtect64



Applications

Dmunicipal drinking water uptakes and distribution
Dhotels and leisure centres
Dhospitals, medical facilities
Dapartment blocks, palaces and public buildings
Dwells

Features

Darray of 64 light barriers to detect fish movements
D30 light barriersto detect immobile fish
Dalarm verification system

Dflow control

- Dtemperature control
- Doverflow detector
- Daccess via LAN, Internet or GSM
- DPCdata-evaluation software
- Ddechlorination option



Benefits

Dcontinuous registration of movements all over the aquarium
Dalarm-relevant signals
Dreduction of false alarms
Dlow maintenance
Deasy service
Dreliable data evaluation
Dadd-ons

Your local bbe dealer...

Envitech Ltd Unit S7, Capital Business Park Parkway. Cardiff CF3 2PU Email: info@envitech.co.uk Tel: +44 2920 364252 Web: www.envitech.co.uk



Specifications

ToxProtect64 (basic instrument without options)

Weight of the instrument	50 kg
Size (H x W x D)	800 x 790 x 444 mm ³
Power supply	110/230V @50/60 Hz
Power consumption	200 W
Sampletemperature	5 - 28 °C
Sampleflow	50 -150 l/h
Turbidity	< 40 FTU
Aquarium volume	91
Housing protection class	IP54
Relayoutputs	4 x 24V1 A contacts
Interface	LAN
Maintenance interval	7 days: 1 h
	7 days, Th
Recommended fish species	zebrafish, tiger barb, minnow, goldfish, bitterling – depending on the temperature; tests for other local fish on request
Recommended fish species	zebrafish, tiger barb, minnow, goldfish, bitterling – depending on the temperature; tests for other local fish on request automatic feeder with selectable intervals
Recommended fish species Feeding of the fish Number of fish	zebrafish, tiger barb, minnow, goldfish, bitterling – depending on the temperature; tests for other local fish on request automatic feeder with selectable intervals 15-20
Recommended fish species Feeding of the fish Number of fish Length of fish	 zebrafish, tiger barb, minnow, goldfish, bitterling – depending on the temperature; tests for other local fish on request automatic feeder with selectable intervals 15-20 4-6 cm
Recommended fish species Feeding of the fish Number of fish Length of fish Recommended tank time of fish	 zebrafish, tiger barb, minnow, goldfish, bitterling – depending on the temperature; tests for other local fish on request automatic feeder with selectable intervals 15-20 4-6 cm 4 weeks
Recommended fish species Feeding of the fish Number of fish Length of fish Recommended tank time of fish Optional features	 zebrafish, tiger barb, minnow, goldfish, bitterling – depending on the temperature; tests for other local fish on request automatic feeder with selectable intervals 15-20 4-6 cm 4 weeks dechlorination system, visual alarm indicator, support stand for desktop operation, pressurereducer/strainer
Recommended fish species Feeding of the fish Number of fish Length of fish Recommended tank time of fish Optional features Optional interfaces	 zebrafish, tiger barb, minnow, goldfish, bitterling – depending on the temperature; tests for other local fish on request automatic feeder with selectable intervals 15-20 4-6 cm 4 weeks dechlorination system, visual alarm indicator, support stand for desktop operation, pressurereducer/strainer 4-20mA (2x), RS232, RS485, GSM, SCADA

biological · biophysical · engineering

bbe moldaenke

Preetzer Chaussee177 | 24222 Schwewntinetal | Germany

 phone
 (+49)
 431
 38040-0

 fax
 (+49)
 431
 38040-10

 e-mail
 info@bbe-moldaenke.de

 net
 www.bbe-moldaenke.de



Part of the EU-funded Techneau project